

KORU Materials List



Device	Component Name	Material(s)
Precision Tubing Sets	Male Cap	Polypropylene (PP)
Precision Tubing Sets	Female Luer Lock	Polyvinyl chloride (PVC), Non-DEHP
Precision Tubing Sets	Tubing	Polyvinyl chloride (PVC), Non-DEHP
Precision Tubing Sets	Female Cap	Polypropylene (PP)
Precision Tubing Sets	Slide Clamp	High-density polyethylene (HDPE)
Precision Tubing Sets	Male Luer	Polyvinyl chloride (PVC), Non-DEHP
Precision Tubing Sets	Bonding Material	Acrylated Urethane
Precision Flow Controller	Tubing	Polyvinyl chloride (PVC), Non-DEHP
Precision Flow Controller	Slide Clamp	High-density polyethylene (HDPE)
Precision Flow Controller	Flow Controller body	"Acrylonitrile Butadiene Styrene (ABS) - Housing Styrene-Ethylene-Butylene-Styrene (SEBS) - Gasket"
Precision Flow Controller	Tubing	Polyvinyl chloride (PVC), Non-DEHP
Precision Flow Controller	Male Luer	Polyvinyl chloride (PVC), Non-DEHP
Precision Flow Controller	Female Cap	Polypropylene (PP)
Precision Flow Controller	Female Luer Lock	Polyvinyl chloride (PVC), Non-DEHP
Precision Flow Controller	Male Cap	Polypropylene (PP)
Precision Flow Controller	Bonding Material	Acrylated Urethane
Precision Flow Controller	Bonding Material	Acrylated Urethane
Precision Flow Controller	Lubricant	Silicone
HlgH-Flo Subcutaneous Needle Sets	Female Luer	Polyvinyl chloride (PVC), Non-DEHP
HlgH-Flo Subcutaneous Needle Sets	Butterfly Cover	Methyl Methacrylate Acrylonitrile Butadiene Styrene (MABS)
HlgH-Flo Subcutaneous Needle Sets	Needle	304 Stainless Steel
HlgH-Flo Subcutaneous Needle Sets	Butterfly Base	Polypropylene (PP)
HlgH-Flo Subcutaneous Needle Sets	Butterfly Wings Closure	Styrene butadiene styrene (SBS)
HlgH-Flo Subcutaneous Needle Sets	Tubing	Polyvinyl chloride (PVC), Non-DEHP
HlgH-Flo Subcutaneous Needle Sets	Needle Protector	Polypropylene (PP)
HlgH-Flo Subcutaneous Needle Sets	Slide Clamp	High-density polyethylene (HDPE)
HlgH-Flo Subcutaneous Needle Sets	Cohesive Tape	APILFCS-10; White kraft paper, Latex Free
HlgH-Flo Subcutaneous Needle Sets	Male Cap	Polypropylene (PP)
HlgH-Flo Subcutaneous Needle Sets	Y-Connector Body	Polyvinyl chloride (PVC), Non-DEHP
HlgH-Flo Subcutaneous Needle Sets	Y-Connector Nut	Polyvinyl chloride (PVC), Non-DEHP
HlgH-Flo Subcutaneous Needle Sets	Adhesive Dressing Tape (OEM product - 3M)	Tegaderm Transparent Film Dressing
HlgH-Flo Subcutaneous Needle Sets	Bonding Material	Acrylated Urethane
Precision Tubing Extension Set	Male Cap	Polypropylene (PP)
Precision Tubing Extension Set	Female Luer Lock	Polyvinyl chloride (PVC), Non-DEHP
Precision Tubing Extension Set	Tubing	Polyvinyl chloride (PVC), Non-DEHP
Precision Tubing Extension Set	Female Cap	Polypropylene (PP)
Precision Tubing Extension Set	Male Luer	Polyvinyl chloride (PVC), Non-DEHP
Precision Tubing Extension Set	Bonding Material	Acrylated Urethane
Y Connector	Y-Connector Body	Polyvinyl chloride (PVC), Non-DEHP
Y Connector	Y-Connector Nut	Polyvinyl chloride (PVC), Non-DEHP