

HIGH-FLO SUBCUTANEOUS SAFETY NEEDLE SETS™

A CUT ABOVE




**High-Flo™
Tricuspid Needle**



**Other brand
Lancet Needle**

LESS TISSUE TRAUMA

Back Cut (Tricuspid) Needle Tip designed to minimize coring and tearing of the skin compared to lancet needles.¹




EVEN FLOW RATE

Needle sets provide even flow to every infusion site.



CUSTOMISABLE

Customize delivery times, number of infusion sites, needle lengths and more.



SAFETY CLOSURE

Butterfly wings snap shut to reduce the risk of accidental needlestick injuries.

LOW RESIDUAL VOLUME (RV) SAVE DRUGS, REDUCE COST

26G Sets Residual Volume

1 needle



RV 0.1 mL

2 needles



RV 0.2 mL

3 needles



RV 0.3 mL

4 needles



RV 0.4 mL

5 needles



RV 0.5 mL

With Y-Connector

7 needles



RV 0.9 mL

8 needles



RV 1.0 mL

24G Sets Residual Volume

1 needle



RV 0.1 mL

2 needles



RV 0.7 mL

3 needles



RV 1.1 mL

4 needles





RV 1.4 mL

AMDEL
MEDICAL



CUSTOMIZED INFUSION SOLUTIONS: CHOOSE THE OPTIMAL NEEDLE LENGTH FOR YOUR PATIENTS' NEEDS



actual 26G needle sizes





26G and 24G

Choose from 4, 6, 9, 12, and 14 mm length options in both 26G and 24G*.

*Not available in 4mm.

ACHIEVE EVEN FLOW RATES WITH UP TO 8 INFUSION SITES SIMULTANEOUSLY USING THE Y-CONNECTOR

Y-connector is sold separately.



Learn more

Reference:

I. Duff C, Ochoa D, Riley P, Murphy E, Zampelli A. Importance of ancillary supplies for subcutaneous immunoglobulin infusion: management of the local infusion site. J Infus Nurs. 2013;36(6):384-390. doi:10.1097/NAN.000000000000009