AMDEL
MEDICAL



FREEDOM60® SYRINGE INFUSION SYSTEM

Reliable infusion therapy, wherever care happens.

A portable, mechanical syringe infusion system for subcutaneous therapy at home and in the community — so treatment fits around real life.

AMDEL MEDICAL × KORU

View the Freedom60

Wherever care happens

The pump travels with the patient, not the ward.



At home

Patients infuse in their own space, at a time that works around work, school and family.



In the community

Suits homecare visits, community clinics and outreach teams delivering therapy off-site.




On the move

No batteries or mains power. Supplied with a travel pouch so treatment travels too.

Mechanical by design

A negator spring drives the Freedom60® — no batteries, no mains power, no electronics to fail. It holds a safe, constant 13.5 psi, and DynEQ™ automatically slows delivery if site pressure rises, so infusions stay steady and predictable.



Why teams choose Freedom60®

Portable and self-contained

A lightweight syringe driver that needs no batteries or mains power. Supplied with a travel pouch so patients keep their independence.

Built for subcutaneous therapy

Used for subcutaneous immunoglobulin and intravenous antibiotics, with volumes from 1ml to 60ml in a standard BD® 50ml syringe.

Safe, constant pressure

Holds an accurate 13.5 psi throughout. DynEQ™ responds automatically by slowing the rate if pressure at the site increases.

Simple to teach and use

Only two controls — wind, load and go. Quick for patients, carers and homecare teams to learn and to support remotely.



Supporting clinicians and homecare teams to deliver infusion therapy beyond the hospital.

Get in touch

Amdel Medical Limited
123a Allerton Road, Mossley Hill, Liverpool L18 2DD

+ 44 (0)151 733 1900
www.amdelmedical.com
Unsubscribe | Privacy Policy | Terms & Conditions