
AMDEL
MEDICAL




FREEDOM60® Syringe Infusion System

This easy-to-use, lightweight, portable pump is for large doses or recurring infusions at home. Used for the infusion of liquid medications, including subcutaneous immunoglobulin and intravenous antibiotics.

AMDEL MEDICAL × KORU

Get in touch

Infuse in 3 Easy Steps



Wind

Wind the large knob back, so the black tab is at the end of its track.



Load

Load prepared syringe.



Go!

Switch ON and go!



FREEDOM60® Administration

Step-by-Step Instructions



Watch on Youtube

FREEDOM60® Key Features

Quality assured

Made in the USA with highest level of quality control for defect-free product. No product maintenance needed.

Safe constant pressure


Highly accurate constant pressure pump for safe delivery. Maintains pressure of 13.5 psi.

DynEQ™ (Dynamic Equilibrium)

If pressure in sites increase, the system immediately responds by automatically slowing down.

Advanced technology

All parts are designed to fit & work together within the pump. The Negator motor is the ‘heart’ of the FREEDOM60®



Stay connected with us to be the first to hear what’s next.

Get in touch

Amdel Medical Limited
123a Allerton Road, Mossley Hill, Liverpool L18 2DD

+ 44 (0)151 733 1900
www.amdelmedical.com
Unsubscribe | Privacy Policy | Terms & Conditions